mpv hotkeys

https://ocv.me/doc/mpv-hotkeys.html

General	Fast seeking	Viewport
Q Quit	Rewind 5 sec	↑ W Magnify (reduce borders)
Q Quit (save position)	→ Forward 5 sec	W Shrink (show full frame)
F Fullscreen	↑ Forward 60 sec	▼ Move viewport
0 Show progress	Rewind 60 sec	▼ Move viewport
S Screenshot		→ Move viewport
S Screenshot (no subtitles)		↑ Move viewport
Open console	Accurate seeking	
ESC Close console	Rewind 1 sec	Zoom out 10%
	Forward 1 sec	Reset zoom
	↑ Forward 5 sec	
Playback	Rewind 5 sec	
Space Play / Pause	Pgdn Previous chapter	Speed
. One frame forward	Pgup Next chapter	[10% slower
, One frame backward] 10% faster
		{ 50% slower
	Subtitles	} 100% faster
Mouse	V Show / Hide	Reset speed
Doubleclick Fullscreen	R Move higher	, total special
Right-click Play/Pause	R Move lower	
Scroll-Up Forward 10 sec	Z Show 0.1 sec earlier	Audio
Scroll-Down Rewind 10 sec	Z Show 0.1 sec later	M Mute
	J Next language	9 Volume Down
	Previous language	O Volume Up
		# Change language

1 of 1 10/29/2020, 2:59 AM